eabxcol

List of tables

- aclgroup
- booking
- container
- container family
- container type
- <u>document</u>
- event
- event type
- identifier_type
- label
- <u>last movement</u>
- metadata attribute
- metadata_schema
- metadata set
- mime type
- movement_type
- multiple type
- object
- object identifier
- object status
- operation
- project
- project_group
- protocol
- <u>sample</u>
- sample metadata
- sample type
- storage
- storage condition
- storage reason
- <u>subsample</u>

aclgroup (Physical Name: aclgroup)

Groupes des logins

Logical Column Name	Physical Column Name	Туре	PK	Nullable
aclgroup_id (PK)	aclgroup_id	INTEGER	PK	NOT NULL
groupe	groupe	VARCHAR(0)		NOT NULL
aclgroup_id_parent	aclgroup_id_parent	INTEGER		

Referenced By

• project group referencing (aclgroup_id)

booking (Physical Name: booking)

Table des réservations d'objets

Logical Column Name	Physical Column Name	Туре	PK	Nullable
booking_id (PK)	booking_id	INTEGER	PK	NOT NULL
uid (<u>FK</u>)	uid	INTEGER		NOT NULL
booking_date	booking_date	TIMESTAMP		NOT NULL
Date de la rése	ervation			
date_from	date_from	TIMESTAMP		NOT NULL
Date-heure de	début de la réservation			
date_to	date_to	TIMESTAMP		NOT NULL
Date-heure de	fin de la réservation			
booking_comment	booking_comment	VARCHAR(0)		
Commentaire				
booking_login	booking_login	VARCHAR(0)		NOT NULL
Compte ayant	réalisé la réservation			
References • object through (uid)			

container (Physical Name: container)

Liste des conteneurs d'échantillon

Logical Column Name	Physical Column Name	Туре	PK	Nullable
container_id (PK)	container_id	INTEGER	PK	NOT NULL
uid (<u>FK</u>)	uid	INTEGER		NOT NULL
container_type_id (FK)	container_type_id	INTEGER		NOT NULL

References

- <u>object</u> through (uid)
- container_type through (container_type_id)

Referenced By

• storage referencing (container_id)

container_family (Physical Name: container_family)

Famille générique des conteneurs

Logical Column Name	Physical Column Name	Туре	PK	Nullable
container_family_id (PK)	container_family_id	INTEGER	PK	NOT NULL
container_family_name	container_family_name	VARCHAR(0)		NOT NULL
is_movable	is_movable	BOOLEAN		NOT NULL

Indique si la famille de conteneurs est déplçable facilement ou non (éprouvette : oui, armoire : non)

Referenced By

• container_type referencing (container_family_id)

container_type (Physical Name: container_type)

Table des types de conteneurs

Logical Column Name	Physical Column Name	Туре	PK	Nullable
container_type_id (PK)	container_type_id	INTEGER	PK	NOT NULL
container_type_name	container_type_name	VARCHAR(0)		NOT NULL
container_family_id (FK)	container_family_id	INTEGER		NOT NULL
storage_condition_id (FK)	storage_condition_id	INTEGER		
label_id (<u>FK</u>)	label_id	INTEGER		
container_type_description	container_type_description	VARCHAR(0)		
Description longue				
storage_product	storage_product	VARCHAR(0)		
Produit utilisé pour	le stockage (formol, alcool))		
clp_classification	clp_classification	VARCHAR(0)		
Classification du ris	que conformément à la direct	tive européenne (CLP	

References

- container family through (container family id)
- storage condition through (storage condition id)
- <u>label</u> through (label id)

Referenced By

- <u>container</u> referencing (container type id)
- <u>sample type</u> referencing (container type id)

document (Physical Name: document)

Documents numériques rattachés à un poisson ou à un événement

Logical Column Nam	Name	Туре		Nullable
document_id (PK)	document_id	INTEGER	PK	NOT NULL
uid (<u>FK</u>)	uid	INTEGER		NOT NULL
mime_type_id (<u>FK</u>)	mime_type_id	INTEGER		NOT NULL
document_import_date	document_import_date	TIMESTAMP		NOT NULL
Date d'import da	ns la base de données			
document_name	document_name	VARCHAR(0)		NOT NULL
Nom d'origine du	document			
document_description	document_description	VARCHAR(0)		
Description libre	du document			
data	data	[-2]		
Contenu du docu	ment			
thumbnail	thumbnail	[-2]		
Vignette au form	at PNG (documents pdf, jpg	g ou png)		
size	size	INTEGER		
Taille du fichier	téléchargé			
document_creation_date	document_creation_date	TIMESTAMP		
Date de création	du document (date de prise	de vue de la photo)	
		-		
References				
• object through (uid)			
	n (mime type id)			

event (Physical Name: event)

Table des événements

Logical Column Name	Physical Column Name	Туре	PK	Nullable
event_id (PK)	event_id	INTEGER	PK	NOT NULL
uid (<u>FK</u>)	uid	INTEGER		NOT NULL
event_date	event_date	TIMESTAMP		NOT NULL

event_type_id (<u>FK</u>)	event_type_id	INTEGER	NOT NULL
still_available	still_available	VARCHAR(0)	
définit ce qu'il	l reste de disponible dar	ns l'objet	
event_comment	event_comment	VARCHAR(0)	
References			

event_type (Physical Name: event_type)

Types d'événement

Logical Column Name	Physical Column Name	Туре	PK	Nullable
event_type_id (PK)	event_type_id	INTEGER	PK	NOT NULL
event_type_name	event_type_name	VARCHAR(0)		NOT NULL
is_sample	is_sample	BOOLEAN		NOT NULL
L'événement s	applique aux échantillons			
is_container	is_container	BOOLEAN		NOT NULL
L'événement s	'applique aux conteneurs			

Referenced By

• event referencing (event_type_id)

identifier_type (Physical Name: identifier_type)

Table des types d'identifiants

Logical Column Name	Physical Column Name	Туре	PK	Nullable
identifier_type_id (PK)	identifier_type_id	INTEGER	PK	NOT NULL
identifier_type_name	identifier_type_name	VARCHAR(0)		NOT NULL
Nom textuel de l	'identifiant			
identifier_type_code	identifier_type_code	VARCHAR(0)		NOT NULL
Code utilisé pour	r la génération des étique	ttes		

Referenced By

• object_identifier referencing (identifier_type_id)

label (Physical Name: label)

Table des modèles d'étiquettes

Logical Column Name	Physical Column Name	Туре	PK	Nullable
label_id (PK)	label_id	INTEGER	PK	NOT NULL
label_name	label_name	VARCHAR(0)		NOT NULL
Nom du moc	dèle			
label_xsl	label_xsl	VARCHAR(0)		NOT NULL
	fichier XSL utilisé pour la graphics.apache.org/fop/)	a transformation FO	Р	
label_fields	label_fields	VARCHAR(0)		NOT NULL
T. 1 1	amps à intégrar dans la Ol	RCODE, séparés pa	r une	

Referenced By

• container_type referencing (label_id)

last_movement (Physical Name: last_movement)

Logical Column Name	Physical Column Name	Туре	PK	Nullable
uid	uid	INTEGER		
storage_id	storage_id	INTEGER		
storage_date	storage_date	TIMESTAMP		
movement_type_id	movement_type_id	INTEGER		
container_id	container_id	INTEGER		
container_uid	container_uid	INTEGER		

metadata_attribute (Physical Name: metadata_attribute)

Table des attributs rattachés à un jeu de métadonnées

Logical Column Name	Physical Column Name	Туре	PK	Nullable
metadata_attribute_id (PK)	metadata_attribute_id	INTEGER	PK	NOT NULL
metadata_set_id (FK)	metadata_set_id	INTEGER		NOT NULL
metadata_schema_id (FK)	metadata_schema_id	INTEGER		
metadata_name	metadata_name	VARCHAR(0)		NOT NULL
Nom de la métad	onnée (creator, name)			
metadata_code	metadata_code	VARCHAR(0)		
Code normalisé d	le la métadonnée (ex : dcte	rms:creator)		
metadata_order	metadata_order	INTEGER		NOT NULL
Ordre d'affichage	des informations dans la g	grille de saisie		
metadata_type	metadata_type	VARCHAR(0)		NOT NULL
metadata_defaultvalue	metadata_defaultvalue	VARCHAR(0)		
metadata_measure_unit	metadata_measure_unit	VARCHAR(0)		
Unité de mesure	utilisée			
metadata_multivalue	metadata_multivalue	BOOLEAN		NOT NULL
metadata_enum	metadata_enum	VARCHAR(0)		

References

- metadata_set through (metadata_set_id)
 metadata_schema through (metadata_schema_id)

Liste des valeurs possibles, séparées par ;

metadata_schema (Physical Name: metadata_schema)

Liste des schémas de métadonnées utilisés

Logical Column Name	Physical Column Name	Туре	PK	Nullable
metadata_schema_id (PK)	metadata_schema_id	INTEGER	PK	NOT NUL
metadata_schema_name	metadata_schema_name	VARCHAR(0)		NOT NUL
Nom complet du schéi	ma			
metadata_schema_short_name	metadata_schema_short_name	VARCHAR(0)		
abréviation habituelle	(CC, DC)			
uri	uri	VARCHAR(0)		
Adresse URI d'accès à	la description du schéma			
Referenced By				

metadata_set (Physical Name: metadata_set)

Jeu de métadonnées permettant de décrire précisément un échantillon

Logical Column Name	Physical Column Name	Туре	PK	Nullable
metadata_set_id (PK)	metadata_set_id	INTEGER	PK	NOT NULL
metadata_set_name	metadata_set_name	VARCHAR(0)		NOT NULL

Referenced By

- <u>sample_type</u> referencing (metadata_set_id)
- metadata attribute referencing (metadata set id)
- <u>sample type</u> referencing (metadata set id)

mime_type (Physical Name: mime_type)

Types mime des fichiers importés

Logical Column Name	Physical Column Name	Туре	PK	Nullable
mime_type_id (PK)	mime_type_id	INTEGER	PK	NOT NULL
extension	extension	VARCHAR(0)		NOT NULL
Extension du f	fichier correspondant			
content_type	content_type	VARCHAR(0)		NOT NULL
type mime off	iciel			

Referenced By

• <u>document</u> referencing (mime_type_id)

movement_type (Physical Name: movement_type)

Type de mouvement

Logical Column Name	Physical Column Name	Туре	PK	Nullable
movement_type_id (PK)	movement_type_id	INTEGER	PK	NOT NULL
movement_type_name	movement_type_name	VARCHAR(0)		NOT NULL
Referenced By				
storage referencingsubsample reference				

multiple_type (Physical Name: multiple_type)

Table des types de contenus multiples

Logical Column Name	Physical Column Name	Туре	PK	Nullable
multiple_type_id (PK)	multiple_type_id	INTEGER	PK	NOT NULL
multiple_type_name	multiple_type_name	VARCHAR(0)		NOT NULL
Referenced By				
• <u>sample_type</u> referencing (multiple_type_id)				

object (Physical Name: object)

Table des objets Contient les identifiants génériques

Logical Column Name	Physical Column Name	Туре	PK	Nullable
uid (PK)	uid	INTEGER	PK	NOT NULL
identifier	identifier	VARCHAR(0)		
Identifiant four	ni le cas échéant par le p	rojet		
object_status_id (FK)	object_status_id	INTEGER		
wgs84_x	wgs84_x	DOUBLE		
Longitude GPS	, en valeur décimale			
wgs84_y	wgs84_y	DOUBLE		
Latitude GPS,	en décimal			
References				

• <u>object_status</u> through (object_status_id)

Referenced By

- <u>sample</u> referencing (uid)
- container referencing (uid)
- event referencing (uid)
- storage referencing (uid)
- booking referencing (uid)
- <u>document</u> referencing (uid)
- object_identifier referencing (uid)

object identifier (Physical Name: object identifier)

Table des identifiants complémentaires normalisés

Logical Column Name	Physical Column Name	Туре	PK	Nullable
object_identifier_id (PK)	object_identifier_id	INTEGER	PK	NOT NULL
uid (<u>FK</u>)	uid	INTEGER		NOT NULL
identifier_type_id (<u>FK</u>)	identifier_type_id	INTEGER		NOT NULL
object_identifier_value	object_identifier_value	VARCHAR(0)		NOT NULL

Valeur de l'identifiant

References

- <u>object</u> through (uid)
- identifier type through (identifier type id)

object_status (Physical Name: object_status)

Table des statuts possibles des objets

Logical Column Name	Physical Column Name	Туре	PK	Nullable
object_status_id (PK)	object_status_id	INTEGER	PK	NOT NULL
object_status_name	object_status_name	VARCHAR(0)		NOT NULL

Referenced By

• <u>object</u> referencing (object_status_id)

operation (Physical Name: operation)

Logical Column Name	Physical Column Name	Туре	PK	Nullable
operation_id (PK)	operation_id	INTEGER	PK	NOT NULL
protocol_id (FK)	protocol_id	INTEGER		NOT NULL
operation_name	operation_name	VARCHAR(0)		NOT NULL
operation_order	operation_order	INTEGER		

Ordre de réalisation de l'opération dans le protocole

References

• protocol through (protocol_id)

Referenced By

• <u>sample_type</u> referencing (operation_id)

project (Physical Name: project)

Table des projets

Logical Column Name	Physical Column Name	Туре	PK	Nullable
project_id (PK)	project_id	INTEGER	PK	NOT NULL
project_name	project_name	VARCHAR(0)		NOT NULL

Referenced By

- <u>sample</u> referencing (project id)
- <u>project group</u> referencing (project_id)

project_group (Physical Name: project_group)

Table des autorisations d'accès à un projet

Logical Column Name	Physical Column Name	Туре	PK Nullable

project_id (PK) (<u>FK</u>)	project_id	INTEGER	PK	NOT NULL
aclgroup_id (PK) (FK)	aclgroup_id	INTEGER	PK	NOT NULL

References

- project through (project_id)aclgroup through (aclgroup_id)

protocol (Physical Name: protocol)

Logical Column Name	Physical Column Name	Туре	PK	Nullable
protocol_id (PK)	protocol_id	INTEGER	PK	NOT NULL
protocol_name	protocol_name	VARCHAR(0)		NOT NULL
protocol_file	protocol_file	BLOB		
Description P	DF du protocole			
protocol_year	protocol_year	SMALLINT		
Année du pro	tocole			
protocol_version	protocol_version	VARCHAR(0)		NOT NULL
Version du pr	otocole			

Referenced By

• operation referencing (protocol_id)

sample (Physical Name: sample)

Table des échantillons

Logical Column Name	Physical Column Name	Туре	PK	Nullable
sample_id (PK)	sample_id	INTEGER	PK	NOT NULL
uid (<u>FK</u>)	uid	INTEGER		NOT NULL
project_id (FK)	project_id	INTEGER		NOT NULL
sample_type_id (<u>FK</u>)	sample_type_id	INTEGER		NOT NULL
sample_creation_date	sample_creation_date	TIMESTAMP		NOT NULL
Date de création	de l'enregistrement dans la	base de données		
sample_date	sample_date	TIMESTAMP		
Date de création	de l'échantillon physique			
parent_sample_id (FK)	parent_sample_id	INTEGER		
multiple_value	multiple_value	DOUBLE		

Nombre initial de sous-échantillons

References

- <u>object</u> through (uid)
- <u>sample</u> through (parent_sample_id)
- project through (project_id)
- <u>sample type</u> through (sample type id)

Referenced By

- <u>sample</u> referencing (sample id)
- <u>sample_metadata</u> referencing (sample_id)
- <u>subsample</u> referencing (sample id)

sample_metadata (Physical Name: sample_metadata)

Logical Column Name	Physical Column Name	Туре	PK	Nullable
sample_id (PK) (FK)	sample_id	INTEGER	PK	NOT NULL
data	data	VARCHAR(0)		NOT NULL
C1				

Champ JSONB pour stockage des données spécifiques de l'échantillon

References

• sample through (sample_id)

sample_type (Physical Name: sample_type)

Types d'échantillons

Logical Column Name	Physical Column Name	Туре	PK	Nullable
sample_type_id (PK)	sample_type_id	INTEGER	PK	NOT NULL
sample_type_name	sample_type_name	VARCHAR(0)		NOT NULL
container_type_id (<u>FK</u>)	container_type_id	INTEGER		
operation_id (<u>FK</u>)	operation_id	INTEGER		
metadata_set_id (<u>FK</u>)	metadata_set_id	INTEGER		
metadata_set_id_second (FK)	metadata_set_id_second	INTEGER		
Second jeu de mé	tadonnées rattaché au type			
multiple_type_id (<u>FK</u>)	multiple_type_id	INTEGER		

multiple_unit	multiple unit	VARCHAR(0)

Unité caractérisant le sous-échantillon

References

- <u>container_type</u> through (container_type_id)
- metadata_set through (metadata_set_id)
- metadata_set_through (metadata_set_id_second)
- operation through (operation_id)
- multiple type through (multiple type id)

Referenced By

• sample referencing (sample_type_id)

storage (Physical Name: storage)

Gestion du stockage des échantillons

Logical Column Name	Physical Column Name	Туре	PK	Nullable
storage_id (PK)	storage_id	INTEGER	PK	NOT NULL
uid (<u>FK</u>)	uid	INTEGER		NOT NULL
container_id (<u>FK</u>)	container_id	INTEGER		
movement_type_id (FK)	movement_type_id	INTEGER		NOT NULL
storage_reason_id (FK)	storage_reason_id	INTEGER		
storage_date	storage_date	TIMESTAMP		NOT NULL
Date/heure d	lu mouvement			
range	range	VARCHAR(0)		
Emplacemen	nt de l'échantillon dans le co	onteneur		
login	login	VARCHAR(0)		NOT NULL
Nom de l'uti	lisateur ayant réalisé l'opéra	ation		
storage_comment	storage_comment	VARCHAR(0)		
Commentair	re			

References

- container through (container id)
- <u>object</u> through (uid)
- movement type through (movement type id)
- <u>storage_reason</u> through (storage_reason_id)

storage_condition (Physical Name: storage_condition)

Condition de stockage

		_		
Logical Column Name	Physical Column Name	Type	PK	Nullable
storage_condition_id (PK)	storage_condition_id	INTEGER	PK	NOT NULL
storage_condition_name	storage_condition_name	VARCHAR(0)		NOT NULL
Referenced By				
• <u>container_type</u> referencing (storage_condition_id)				

storage_reason (Physical Name: storage_reason)

Table des raisons de stockage/déstockage

Logical Column Name	Physical Column Name	Туре	PK	Nullable		
storage_reason_id (PK)	storage_reason_id	INTEGER	PK	NOT NULL		
storage_reason_name	storage_reason_name	VARCHAR(0)		NOT NULL		
Referenced By • storage referencing (storage reason id)						
<u> </u>	, (000111130_10110111_101)					

subsample (Physical Name: subsample)

Table des prélèvements et restitutions de sous-échantillons

Logical Column Name	Physical Column Name	Туре	PK	Nullable
subsample_id (PK)	subsample_id	INTEGER	PK	NOT NULL
sample_id (<u>FK</u>)	sample_id	INTEGER		NOT NULL
subsample_date	subsample_date	TIMESTAMP		NOT NULL
Date/heure de 1	'opération			
movement_type_id (FK)	movement_type_id	INTEGER		NOT NULL
subsample_quantity	subsample_quantity	DOUBLE		

Quantité prélev	ée ou restituée		
subsample_comment	subsample_comment	VARCHAR(0)	
subsample login	subsample login	VARCHAR(0)	NOT NULL

Login de l'utilisateur ayant réalisé l'opération

References

- sample through (sample_id)movement_type through (movement_type_id)